



tierratech[®]
ULTRASONIC CLEANING SYSTEMS

LT-150 PRO

(5,3 LITERS)



Capacidad: 5,3 liters

Internal dimensions: 295x145x152 mm

External dimensions: 328x177x318mm

Weight: 6 Kg

ELECTRICAL CHARACTERISTICS

- Power Supply: 220V.
- Ultrasonic power: 150W (300W p-p).
- Heating resistance: 150 W.

ULTRASONIC CHARACTERISTICS

- Ultrasonic power: 150W (300W p-p).
- Working frequency: 40 KHz. Frequency scanning system (sweep system $\pm 2\%$).
- 3 high-performance IBL piezo-electric transducer on steel-titanium.
- Digital control board:
 - Main switch.
 - Ultrasonic time selection: from 0 to 60 min.
 - Temperature selection: from 0° to 80°C.

GENERAL CHARACTERISTICS AND COMPLEMENTS

- Tank made in stainless steel AISI 316.
- External panel made in print-free steel.
- Ball valve (1/2") for draining the water.
- 1 stainless steel basket.
- Stainless steel cover.

GUARANTEE AND QUALITY:

Every equipment manufactured by TIERRA TECH in conformed European Directives about CE certificate. Our machines are due to rules about electric and magnetic compatibility and security of electrical components: Directives: 2004/108/CEE, 2006/95/CEE, EN 55011, EN 60204, EN 61010. This unit is guaranteed for 1 year in case of any manufacture fault.